



ELECTRIC HEATING INCUBATOR



BXP-N series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

ELECTRIC HEATING INCUBATOR >>

BXP-N series

Product Features >

- LCD display, real-time display of working status;
- BRIGHT I control system, automatically compensate parameters based on the environment;
- With countdown reservation function, the reservation time range is 99 hours 59 minutes;
- Programmable design: 10 segments 100 cycles program setting;
- Double over temperature alarm function (A. controlled by temperature sensor system, B. independent temperature control switch);
- The 6-stage speed regulation of the fan,(BXP-N series no fan);
- Preheating chamber design, to ensure rapid heating and good heat distribution effect;
- Double door design, with glass inner door;
- Equipped with 1 RS-485 interface.



Lock



Fan



Shelf



| Model | BXP-70N | BXP-140N | BXP-300N |
|------------------------------------|-------------|--------------|---------------|
| Volume [L] | 70 | 140 | 300 |
| Temp. Range [°C] | | RT+5~100 | |
| Temp. Fluctuation at 37°C [°C] | | ±0.5 | |
| Temp. Uniformity at 37°C [°C] | | ±0.6 | |
| Power [W] | 300 | 480 | 800 |
| Internal Dimensions W×D×H [mm] | 410×390×450 | 570×440×570 | 720×570×745 |
| Housing Dimensions W×D×H [mm] | 575×566×776 | 735×616×896 | 885×746×1071 |
| Number of Shelves (std./max.) [pc] | 2/5 | 2/7 | 3/9 |
| Net Weight /Gross Weight [kg] | 46/64 | 57/75 | 97/118 |
| Package Size [mm] | 695*726*926 | 855*776*1046 | 1005*906*1221 |